ISTRUZIONI D'USO TWINFUSER

Destinazione d'uso

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INSTRUCTION FOR USE TWINFUSER

Intended use

Device intended to exert pressure on infusion bags in the medical field both with adult and pediatric patients.

The use of two bags at once allows to perform treatment on two patients concurrently, to use two different bags on the same patient or to always have a spare bag ready to be used in case of emergency.

Twinfuser Principle of Operation:

Twinfuser pressure infusion bag is made of an internal expandable bag (also known as inflatable cushion) between the infusion bags. There is one wraparound band which at the same time wrap fluid bags and cushion.

The entire procedure of bag loading and system preparation can be easily done by one operator over an horizontal surface, which cannot be done with other devices made by competitors because they require the presence and intervention by two operators at once.

It is a completely autonomous device which does not require connections to electrical power supplies. It just requires a support made by an I.V. of any kind which can be commonly found into every healthcare facility.

Contraindications

None if used for the declared purpose.

Instructions for use

TWINFUSER T500

To be used with 1 500cc. bag

To be used with 1 500cc. bag + 1 500cc. bag (tandem configuration)

TWINFUSER T1000

To be used with 1 1000cc. bag

To be used with 1 1000cc. bag + 1 1000cc. bag (tandem configuration)

To be used with 1 1000cc. bag + 1 500cc. bag (tandem configuration)

TWINFUSER T3000

To be used with 1 3000cc. bag

To be used with 1 3000cc. bag + 1 3000cc. bag (tandem configuration).

Position Twinfuser on a flat surface with the tube and feed pump facing the operator.

Open the zip fastener thus making the two compartments to house the infusion fluid bags available. Make sure that the tapes (7) pass through the central eyelets of each bag.

Close the system with the zip fastener. The system can now be hung on a stand.

Check that the screw knob (1) is closed.

Open the 3-way cock (2) and check that the bag (5) is in communication with the pump (8) as shown in Fig. 1. Pressurise Twinfuser with the pump (8) until reaching the desired pressure.

Close the clamp (3) to prevent undesired system leaks and close the 3-way cock (2) as shown in Fig. 2.

The pressure gauge shows an approximate value, expressed in mmHg, of the internal pressure of the inflatable cushion. If a pressure of 300 mmHg is exceeded, the pressure gauge will indicate a red alarm area and automatically activate the overpressure valve.

To quickly depressurise the system, open the clamp (3) and the 3-way cock (2) as shown in Fig. 3.

Warnings

If present on the label. Caution: this product contains phthalate. There is evidence suggesting that phthalate exposure from medical procedures may be harmful for the treatment of children, pregnant women or nursing women. Use of this product does not represent a high exposure risk to phthalate, however, as a precautionary measure its use should be limited to essential procedures.

The device must be used by qualified medical/nursing staff.

Non-sterile. Do not sterilize.

Reusable.

During the use of the device wear gloves.

The device may be reused on the condition that it is appropriately cleaned before each reuse. It can be cleaned by rubbing with normal neutral detergents.

Maximum number of uses: 10 times.

Expiry: 5 years provided that the packaging is undamaged and if stored in normal storage conditions (-20°/+50°C).

They are potentially contaminated wastes which shall be disposed of according to law in force as regards healthcare facilities.

The pressure measurement indicated by the device is approximate. Do not use this system to check the pressure values of the patient.

When inflating the system always check on the pressure gauge connected to the device that the safety pressures are not exceeded.

During the infusion procedure keep monitored the pressure indicated by the manometer and, in case it has been decreased, inflate the pump up to maximum 300 mmHg.

In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.

The choice of inflate pressure and deflation quickness, which depend upon the kind of surgery treatment and the used fluid, are under responsibility of the user. Do not attempt to open the zip fastener when Twinfuser is under pressure.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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