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Reusable non-invasive ventilation (NIV) mask

INTENDED USE

Reusable mask for CPAP or Bilevel ventilatory support.

INDICATIONS FOR USE

Device indicated for the treatment of respiratory insufficiency treatable with CPAP or Bilevel non-invasive ventilation methods.

The mask is intended for reuse on the same patient or on different patients after reconditioning in the manner described below.

LIMITS OF USE

Use in hospital or nursing facility settings only.

The device must be used by trained and qualified medical/nursing staff.

Use with monitoring lines suitable for connection to 5mm pressure ports.

The cap for insertion of nasogastric tubes (SNG cap), if present, can be used for tubes with a diameter of ≥ 12 Ch.

The mask is intended for re-use on several patients according to the disinfection or sterilisation procedure:

TYPE OF PROCEDURE	PATIENT TYPOLOGY
High-level manual disinfection with CIDEX OPA	adult and paediatric patients
High-level manual disinfection ANIOXIDE 1000	
Sterilisation with hydrogen peroxide plasma gas	
Ethylene Oxide (ETO) sterilisation	adult and paediatric patients weighing > 30 kg

The mask is intended for reuse on several patients, both adults and children weighing >30 kg.

Do not use for the administration of anaesthetic gases or medicinal products.

NOTE: for use of the mask with only angled connectors without anti-asphyxiation valve:

Use only with ventilators equipped with adequate alarms and safety systems able to detect interruptions in ventilation. This mask configuration does not include anti-asphyxiation safety systems that allow the patient to breathe even in case of ventilator failure.

CONTRAINDICATIONS

- The device must not be used on patients who are uncooperative or unconscious and cannot remove the mask on their own.
- The device is not suitable for use on patients taking medications that may cause vomiting, patients with nocturnal vomiting, hiatal hernia or functional damage to the cardiac sphincter.
- This mask is not suitable for use in patients with the following conditions: glaucoma, recent eye surgery or eye dryness.

COMPLICATIONS/SIDE EFFECTS

Shortness of breath

Gastric dilatation

Eye infections

CO₂ retention (rebreathing)

Facial trauma

Skin irritation

Eye and airway dryness

Pneumothorax

Claustrophobia

The following must be reported to the treating physician:

unusual chest discomfort, eruption, severe headache, eye dryness or pain, blurred vision while using the device.

WARNINGS AND PRECAUTIONS

In case of an evident change in performance or inadequate performance of the device, it is recommended to replace it.

Monitor the patient's vital signs during therapy.

If there is evidence of worsening of the patient's condition or no improvement within the prescribed time frame, alternative ventilatory techniques should be evaluated.

When used with humidified ventilation gases, check that any condensation does not cause discomfort to the patient.

It is the responsibility of the medical/nursing staff to select the device most suitable for the patient in terms of configuration and size; the appropriate choice will ensure good seal.

Use on intact skin.

Expiry: 5 years in intact packaging and when stored in normal storage conditions.

Non-invasive ventilation should be discontinued in the presence of:

- worsening of the state of consciousness and respiratory distress
- loss of airway protection
- PaCO₂ remaining unchanged (in 2 successive EGAs performed at a maximum distance of 1 hour)
- persistent severe hypoxaemia
- severe and uncontrollable haemodynamic instability
- patient/ventilator synchronisation problems
- uncontrollable secretions
- interface intolerance

Mask with standard unventilated angled connector (BLUE) or with unventilated angled connector with bronchoscopy access (GREEN)

This mask configuration is not equipped with an integrated opening for exhaled gas escape, therefore, the breathing circuit needs to have an adequate expiratory device. Use the breathing circuit recommended by the manufacturer of the ventilation equipment.

Use only with a ventilator equipped with adequate alarms and safety systems against possible ventilator malfunction. This mask configuration does not include an anti-asphyxiation valve that would allow the patient to breathe in the event of a ventilator malfunction.

- This mask is not suitable to provide assisted ventilation as vital support.
- To minimize the risk of vomiting during sleep, avoid eating or drinking for three (3) hours before use.
- Check that the mask and harness are the right size for the patient. In case of excessive discomfort, reposition the mask and the harness.

Mask with ventilated angled connector with anti-asphyxiation valve (ORANGE)

This mask configuration is equipped with an integrated opening for exhaled gas escape; this opening represents continuous leakage from the breathing circuit.

Do not obstruct!

- This mask is not suitable to provide assisted ventilation as vital support.
- To minimize the risk of vomiting during sleep, avoid eating or drinking for three (3) hours before use.
- Check that the mask and harness are the right size for the patient. In case of excessive discomfort, reposition the mask and the harness.
- A minimum pressure of 3.0 cmH₂O is required.
- This mask should only be used with positive airway pressure devices according to the indications of the treating physician. Only use the mask when the positive airway pressure device is turned on and working properly. The integrated openings for exhaled gas escape must never be obstructed.

INSTRUCTIONS FOR USE

Inspection and functional test before use

- Before first use and before every reuse, check that there is no damage or deterioration that may affect proper functioning of the device.
 - Check the angled connectors for deformation, breakage, leakage and obstruction.
 - Check that the angled connector is attached to the mask and that it is properly connected to the breathing circuits.
 - Clean the patient's face thoroughly before use.
 - Select the appropriate mask size; the mask seal must be fitted so that it sits above the eyebrows and below the lower lip. (Figure 1). Keep the harness in front of the mask and place the lower part of the seal below the patient's lower lip. Pull the harness over the patient's head and place the mask on the patient's face. Make sure that the seal rests below the patient's lower lip and above the eyebrows. Adjust the nuchal on the top of the patient's head and forehead while holding the mask in place. Attach both the lower harness clips to the mask; this should produce an audible click. Adjust both the upper harness straps evenly to fit the patient's head. Once the upper straps are fixed, adjust the lower straps evenly. Do not over-tighten the harness. Follow these steps again each time the mask is removed or repositioned.
- To remove the mask from the patient, detach the lower harness clips from the mask, move the harness from the head and lift the mask from the face.


Figure 1

For models equipped with anti-asphyxiation valve

- Check functioning of the anti-asphyxiation valve on the angled connector.
- Check that when the ventilator is off, the anti-asphyxiation valve membrane is positioned so as to allow the patient to breathe ambient air from the angled connector openings.
- Turn on the ventilator; the anti-asphyxiation valve membrane must close so that the patient receives the flow from the ventilator.
- If the valve does not function as indicated, replace the angled connector.
- Do not obstruct the anti-asphyxiation valve opening.
- Always check the anti-asphyxiation valve before use for any malfunction or obstruction caused by the patient's secretions.

Use of the interchangeable angled connectors


Remove the mask from the patient.

Remove the existing angled connector from the mask by holding the ring firmly on the mask and turning the angled connector counterclockwise to disconnect the bayonet coupling.

Remove the angled connector.

Place a new angled connector in the ring on the mask and activate the bayonet coupling by turning it clockwise until it fits into place.

Do not force or use tools.

Use of the unventilated angled connector with bronchoscopy access (GREEN)

The unventilated angled connector with bronchoscopy access is an accessory to provide the patient with CPAP and Bilevel therapy during a bronchoscopy procedure.

This accessory should only be used temporarily during bronchoscopy procedures.

WARNINGS ON REUSE

The reconditioning instructions in this leaflet must be followed before each reuse of the device on the same or a new patient.

Do not use gamma rays or ultrasonic cleaning equipment. Do not sterilize with dry or wet steam sterilizers. Do not sterilize or disinfect with components or derivatives of phenol and alkylamine (glucoprotamine).

Do not clean with solvent- or chloride-based solutions.

Do not use hard brushes or other materials that could damage the silicone parts of the device. Avoid contact with hard or sharp objects.

The cleaning operations must always be carried out prior to the next steps of disinfection or sterilization.

After cleaning and subsequent disinfection or sterilization, visually inspect the product. Slight yellowing/discoloration is considered normal.

Do not reuse if components are clearly cracked or separated.

If the cleaning procedure does not produce adequate results, disinfection or sterilization of the product is not recommended, but it should be disposed of as described in the dedicated section.

Properly store the mask after disinfection or sterilization to avoid contamination.

Refer to the technical data sheets of the disinfectants for the relative operating procedures and safe handling.

The harness, angled connector, SNG cap (if present) and monitoring port caps must be separated from the mask and from the angled connector before cleaning.

The harness and SNG cap are single patient components and must therefore be replaced before each reuse on the same or a new patient.

LIMITS OF REUSE

Device sold in a non-sterile configuration. Reconditionable and reusable up to a maximum of 10 times.

After a cleaning phase, the device may be reprocessed by high-level disinfection (CIDEX OPA for mask and connectors, ANIOXIDE 1000 for mask only) or sterilization (ETO or hydrogen peroxide gas plasma).

RECONDITIONING INSTRUCTIONS

Point of use

See the cleaning instructions below.

Content and transport

It is recommended to recondition all the products as soon as possible after each use and before subsequent reuse.

Preparation for cleaning

Prior to cleaning and subsequent disinfection or sterilization, separate the angular connector, SNG cap (if present), harness and monitoring port caps from the mask.

The harness and SNG cap cannot be treated with the disinfection or sterilization methods described, as these are single patient accessories that must always be replaced before each reuse.

Cleaning is required before each reconditioning. Before disinfection or sterilization, clean the mask and connectors following the instructions below (manual cleaning).

Manual cleaning:

Wash your hands with a mild detergent, such as a neutral pH liquid soap or trisodium phosphate. Rinse your hands thoroughly with clean water and dry with a soft absorbent material that does not release fibres. Dismantle the mask by removing the harness, angled connector, SNG cap (if present) and monitoring port caps.

Immerse the mask, angled connector, SNG cap (if present) and monitoring port caps in commercially available anionic detergent for 5 minutes, for example DIALZIMA ULTRA diluted in water (0.15%).

Rinse for 10 seconds under running water.

Do not use hard brushes or other materials that could damage the products. Do not use washing machines. Dry at room temperature.

Do not dry on a hot surface. Air dry avoiding direct exposure to sunlight.

WARNING

The complete procedure includes the steps "preparation for cleaning", "manual cleaning" and subsequent disinfection or sterilization.

After cleaning, proceed with disinfection or sterilization prior to subsequent reuse on the same or a new patient.

High-level manual disinfection

Use one of the disinfecting agents listed below. Respect the concentrations and immersion times indicated by the manufacturer.

Rinse thoroughly with clean water and dry with a soft absorbent material that does not leave fibres on the devices. Dry at room temperature.

Treatment	Exposure time	Active ingredient	Solution	Rinsing	Maximum disinfection cycles
CIDEX OPA	5 minutes in immersion	0.55% ortho-phthalaldehyde	/	10 minutes immersion in water for 3 times	10
ANIOXIDE 1000 (only for mask and not for connectors)	5 minutes in immersion	0.15% peracetic acid	Activation as per the manufacturer's data sheet	1 minute immersion in water for 1 time	10

Sterilization:

Intersurgical SpA has verified the effectiveness of sterilization on the medical devices covered by these Instructions for Use. The parameters given below are to be considered recommended. Any different cycles must be validated by the user.

After sterilization, respect the degassing times (for ETO) and finally let the devices dry at room temperature.

Sterilization with hydrogen peroxide gas plasma

Treatment	Exposure time	Sterilization temperature	Sterilizing agent	Minimum value of vacuum	Maximum disinfection cycles
STERRAD	See machine manual	45°C-55°C	Hydrogen Peroxide >58%	< 80 Pa	10

Sterilization with ethylene oxide (ETO)

Treatment	Exposure time	Range temperature during gas exposure time	Relative Humidity	Gas composition	Maximum disinfection cycles
Ethylene Oxide	18 hours + 6 days degassing	45°C (± 5°C)	52%	15% ETO and 85% CO ₂	10

WARNING: USE THE MEDICAL DEVICE NOT EARLIER THAN 8 DAYS AFTER STERILIZATION.

STORAGE

The device must be stored at a temperature between -20° and +50°C.

The manufacturer of this device validates and provides cleaning, disinfection and sterilization instructions in its capacity as manufacturer of reusable medical devices. The reconditioner (the user who carries out the cleaning, disinfection and sterilization procedures) is responsible for ensuring that reconditioning is carried out in a dedicated facility with equipment, materials and personnel suited to achieve the desired result. This normally requires validation and continuous monitoring of the reconditioning process. Similarly, any deviation by the reconditioner from the instructions provided by the manufacturer should appropriately be evaluated by the reconditioner in terms of effectiveness and potential adverse consequences.

OPERATING SPECIFICATIONS

+5° to +40°C in a humidity range of 15-95%

DEVICE DISPOSAL

Dispose of the devices as potentially hazardous material in accordance with hospital legislation.

IT IS RECOMMENDED TO KEEP THESE INSTRUCTIONS.

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority where the user and/or patient are located.

LEGEND OF SYMBOLS

MOD	Model/size
	Code
	Lot
	Expiry date
QTY	Quantity
	See instructions for use attached to the device
	Warning
	Manufacturer
	Do not use if the package is damaged
	Non-sterile
	Fragile
	Temperature limitations
	Do not open the packaging with a knife
	The CE marking includes the TUV Rheinland Italia (Notified Body) identification number. Product compliant with MDR (EU) 2017/745 requirements
	Keep away from sunlight
	Keep away from rain
	Phthalate-free
	Latex-free
	Medical device

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