(it) ACCESSORIO PER MASCHERE

(en) ACCESSORY FOR MASKS

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Istruzioni per l'uso / Instructions for use / Mode d'emploi / Gebrauchsanleitung

(it)

RENDIMENTI

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(en)

PERFORMANCE

REQUIREMENT		SETTING	TO PATIENT
0 % (± 5%)	PEEP cmH ₂ O	O ₂ FLOW Ⓐ + Ⓑ	FLOW L/min (± 8%)
100%	For 100% oxygen treatment connect both tubes to oxygen sources and close the Venturi valve. Feed with a minimum oxygen flow of 40 L/min		

Destinazione d'uso

The accessory is used in combination with masks for CPAP (Continuous Positive Airway Pressure) treatment. To be used in hospital environment when high flow of air/oxygen mixture at variable and controlled O₂% from 30% to 100% is required.

Indications for use

Indicated for the treatment of acute respiratory insufficiency of hypoxemic patients, and post-surgery atelectasis.

Limits of use

The device must be used by qualified and trained medical/nursing staff. To be used for administration of air and oxygen. To be used on adult patients.

If used properly, the device can be used continuously for a maximum of 7 days after which it must be replaced.

To be used solely with oxygen flowmeters capable of delivering the flows indicated in the table. With one O₂ source a maximum FiO₂ of 30% is obtained; therefore, to obtain a greater FiO₂ use two O₂ sources.

Use flowmeters with a 7mm diameter hose tail.

To be connected with respiratory systems/patient interfaces featuring 22mm connections (compliant with ISO 5356-1 standards).

Where it can be used

The accessory allows the delivery of high flow CPAP therapy using compressed oxygen sources (wall-mounted distribution systems or cylinders) available in hospital departments and ambulances.

HOW TO USE ACCESSORY

1. Remove the accessory from the package. Connect the device to the mask/patient interface gas inlet connector. Verify safe connections.
2. Attach the oxygen input connector/s A and B to the flowmeter/s (see figure 1). Check proper connection before activating the gas supply.


Fig. 1

3. Activate the flowmeter/s according to the performance table attached to the device.
4. To supply O₂% equal to 100% the Venturi valve has to be inactivated by sliding the transparent cover (see figure 2a) until the air inlet is occluded (lock to closed position – see figure 2b).
5. Only when the air inlet/transparent cover is closed, the total flow to the patient is equal to the sum of the oxygen flows from the two flowmeters. To restore the operation of the Venturi valve, slide and lock the transparent cover to the opened position (figure 2a) and adjust therapy conditions according to the performance table.

NOTE: The cover must be fully engaged to the opened/closed positions.


Fig. 2a

Fig. 2b

The manometer provides an indication of the presence of pressure within the breathing system, the indication is however approximate.

WARNINGS AND CAUTIONS

Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, the materials of the device could be damaged by disinfection methods, the product performances may be impaired developing a safety risk for patients.

1. Depending on the clinical condition of the patient the medical/nursing staffs are responsible for defining the need for more frequent replacement of the device.
2. It is recommended to use the settings provided by Intersurgical within the performance table attached to the device. These settings ensure sufficient flows for proper CPAP therapy. Higher flow settings would generate unnecessary oxygen consumption and excessive noise.
3. Always check that the oxygen supply lines are safely connected. In the event of oxygen tubes disconnection from the device, patient ventilation will fail and oxygen could be released into the environment.
4. In the events of an evident change in performance or inadequate performance of the device, it is recommended to replace it.
5. It is recommended to monitor the clinical parameters of the patient during the therapy.
6. If the patient's conditions deteriorate or there is no improvement within the expected time, it is recommended to evaluate alternative ventilation techniques.
7. Do not use oxygen bubble humidifiers with the gas supply flowmeters.
8. Non-sterile. Do not sterilize.
9. Expiry: 5 years provided that the packaging is undamaged and the product is stored in normal storage conditions (-20°C / +50°C).
10. Dispose of the device immediately after use in compliance with the current laws and regulations.
11. Use only with devices that features connections in compliance with the ISO 5356-1 standards.
12. Caution: the device contains metal; do not use in departments where the presence of metal is a source of risk to the safety of the patient or third parties.
13. The manometer provides an indication of the presence of pressure within the breathing system, the indication is however approximate.
14. Do not use on another patient.

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority where the user and/or patient are located.

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