
VENTUPLUS QUIETO

- (it) ACCESSORIO PER CASCHI E MASCHERE**
- (en) ACCESSORY FOR HELMETS AND MASKS**
- (fr) ACCESSOIRES POUR CASQUES ET MASQUES**
- (de) ZUBEHÖR FÜR HELME UND MASKEN**
- (es) ACCESORIO PARA CASCOS Y MASCARILLAS**
- (pt) ACESSÓRIOS PARA CAPACETES E MÁSCARAS**
- (no) TILBEHØR FOR HJELMER OG MASKER**
- (fi) LISÄLAITE KYPÄRIIN JA MASKEIHIN**
- (sv) TILLBEHÖR TILL HJÄLMAR OCH MASKER**
- (da) TILBEHØR TIL HJELME OG MASKER**
- (el) ΒΟΗΘΗΤΙΚΟ ΕΞΑΡΤΗΜΑ ΓΙΑ ΚΑΣΚΕΣ ΚΑΙ ΜΑΣΚΕΣ ΑΕΡΙΣΜΟΥ**
- (pl) AKCESORIA DO KASKÓW I MASEK**
- (sl) PRIPOMOČEK ZA ČELADE IN MASKE**
- (et) KIIVRITE JA MASKIDE TARVIK**
- (ro) ACCESORII PENTRU CĂȘTI ȘI MĂȘTI**
- (hr) DODATAK ZA KACIGE I MASKE**

*(it) Istruzioni per l'uso - (en) Instructions for use - (fr) Instructions d'utilisation - (de) Gebrauchsanweisung -
(es) Instrucciones de uso - (pt) Instruções de utilização - (no) Instruksjoner for bruk - (fi) Käyttöohjeet - (sv) Bruksanvisning -
(da) Brugsanvisning - (el) Οδηγίες χρήσης - (pl) Instrukcja stosowania - (sl) Navodila za uporabo - (et) Kasutusjuhend -
(ro) Instrucțiuni de utilizare - (hr) Upute za upotrebu*

(it) NOTA: suggerimento di configurazione - (en) NOTE: suggested configuration - (fr) REMARQUE : configuration suggérée - (de) HINWEIS: vorgeschlagene Konfiguration - (es) NOTA: Configuración sugerida - (pt) NOTA: configuração sugerida - (no) OBS: foreslått konfigurasjon - (fi) HUOMAUTUS: ehdotettu kokoonpano - (sv) ANMÄRKNING: Förslag på konfiguration - (da) BEMÆRK: foreslået konfiguration - (el) ΣΗΜΕΙΩΣΗ: συνιστώμενος τύπος - (pl) UWAGA: konfiguracja sugerowana - (sl) OPOMBA: predlagana konfiguracija - (et) MÄRKUS: Soovitatud seadistus - (ro) NOTĂ: configurație recomandată - (hr) NAPOMENA: predložena konfiguracija

FIG. 1 / ABB. 1 / KUVÄ 1 / ILL. 1 / EIK. 1 / RYS. 1 / SL. 1 / JN 1 / SL. 1


FIG. 2 / ABB. 2 / KUVÄ 2 / ILL. 2 / EIK. 2 / RYS. 2 / SL. 2 / JN 2 / SL. 2


FIG. 3 / ABB. 3 / KUVÄ 3 / ILL. 3 / EIK. 3 / RYS. 3 / SL. 3 / JN 3 / SL. 3



FIG. 4 / ABB. 4 / KUVÄ 4 / ILL. 4 / EIK. 4 / RYS. 4 / SL. 4 / JN 4 / SL. 4

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(en) CONFIGURATIONS:

VENTUPLUS QUIETO CP932 WITH DUAL OXYGEN FLOWMETER, SUGGESTED CONFIGURATION FOR O2% ABOVE 40% (UP TO 100%) SEE FIG. 1

VENTUPLUS QUIETO CP932 WITH TWO SINGLE OXYGEN FLOWMETERS, SUGGESTED CONFIGURATION FOR O2% ABOVE 40% (UP TO 100%) SEE FIG. 2

VENTUPLUS QUIETO CP932 WITH SINGLE OXYGEN FLOWMETER, SUGGESTED CONFIGURATION FOR O2% EQUAL TO 40% SEE FIG. 3

VENTUPLUS QUIETO CP934 WITH DUAL OXYGEN FLOWMETER, SUGGESTED CONFIGURATION FOR O2% ABOVE 40% (UP TO 70%) SEE FIG. 4

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

CP932 PERFORMANCE

REQUIREMENT		SETTING	TO PATIENT
O ₂ (± 5%)	PEEP cmH ₂ O	O ₂ FLOW	FLOW L/min (± 8%)
100%	For 100% oxygen treatment connect both tubes to oxygen sources and close the Venturi valve. Feed with a minimum oxygen flow of 40 L/min		

CP934 PERFORMANCE

REQUIREMENT		SETTING	TO PATIENT
O ₂ (± 5%)	PEEP cmH ₂ O	O ₂ FLOW	FLOW L/min (± 8%)

INTRODUCTION

The VENTUPLUS accessory is based on the Venturi tube operating principle: a constant oxygen flow blows through the tube which has a constriction (nozzle). The speed and pressure of the oxygen flow undergo a variation which allows aspirating environmental air from the area surrounding the Venturi tube. The total flow is a mixture of air and oxygen whose O₂% can be modified by using the additional oxygen source if present. The performance table attached to the product indicates how to obtain the combinations of O₂% and total patient flow at each PEEP level.

DESTINATION OF USE

The accessory is used in combination with helmets and masks for CPAP (Continuous Positive Airway Pressure) treatment. To be used in hospital environment when high flow of air/oxygen mixture at variable O₂% is required.

INDICATIONS FOR USE

To be used with helmets and masks for CPAP (Continuous Positive Airway Pressure) treatment. It can deliver a high flow of air/oxygen mixture according to the setting by healthcare personnel.

RESTRICTIONS OF USE

- The device must be used by qualified and trained medical/ nursing staff.
- To be used for administration of air and oxygen.
- If used properly, the device can be used continuously for a maximum of 7 days after which it must be replaced.
- To be used solely with oxygen flowmeters capable of delivering the flows indicated in the table
- Use only with oxygen flowmeters featuring 7mm diameter hose tail.
- To be connected with respiratory systems/patient interfaces featuring 22mm connections (compliant with ISO standards).

WHERE IT CAN BE USED

The VENTUPLUS accessory allows the delivery of high flow CPAP therapy using compressed oxygen sources (wall-mounted distribution systems or cylinders) available in hospital departments and ambulances.

HOW TO USE VENTUPLUS QUIETO

1. Remove the VENTUPLUS QUIETO accessory from the package. Connect the device to the breathing system/patient interface gas inlet connector. Verify safe connections.
2. Attach the oxygen input connector/s A and B to the flowmeter/s (see FIGURE 5). Check proper connection before activating the gas supply.


Fig. 5

3. Activate the flowmeter/s according to the performance table attached to the device.
4. To supply O₂% equal to 100%, apply the cap present in the package to the air inlet filter (see FIGURE 6).


Fig. 6

13. The efficacy of the therapy is considerably affected by the PEEP value set. Therefore, the medical staff are advised to carefully evaluate the pressure level most suited to the clinical condition of the patient. Using too low pressure may not be sufficient for alveolar recruitment. Using too high pressure may be a source of alveolar over distension.



Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

WARNINGS AND CAUTIONS

⊗ Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, the materials of the device could be damaged by disinfection methods, the product performances may be impaired developing a safety risk for patients.

1. Depending on the clinical condition of the patient the medical/nursing staffs are responsible for defining the need for more frequent replacement of the device.
2. It is recommended to use the settings provided by Intersurgical within the performance table attached to the device. These settings ensure sufficient flows for proper CPAP therapy. Higher flow settings would generate unnecessary oxygen consumption and excessive noise.
3. Always check that the oxygen supply lines are safely connected. In the event of oxygen tubes disconnection from the device, patient ventilation will fail and oxygen could be released into the environment.
4. In the events of an evident change in performance or inadequate performance of the device, it is recommended to replace it.
5. It is recommended to monitor the clinical parameters of the patient during the therapy.
6. If the patient's conditions deteriorate or there is no improvement within the expected time, it is recommended to evaluate alternative ventilation techniques.
7. Do not use oxygen bubble humidifiers with the gas supply flowmeters.
8. Non-sterile. Do not sterilize.
9. Expiry: 5 years provided that the packaging is undamaged and the product is stored in normal storage conditions (-20°C / +50°C).
10. Dispose of the device immediately after use in compliance with the current laws and regulations.
11. Use only with devices that features connections in compliance with the ISO 5356-1 standard.
12. The pressure relief valve (if present) allows not to exceed 25 cm H₂O of the helmet internal pressure.

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Intersurgical SpA - Via Morandi, 12 - 41037 Mirandola (MO) Italy