



PEP Therapy Range Positive Expiratory Pressure




Oxygen and Aerosol Therapy ▪ Respiratory Physiotherapy



The Intersurgical PEP therapy range

Intersurgical offer a range of products for the delivery of positive expiratory pressure therapy. The modular design of the product range enables clinicians to select from nine colour coded PEP resistors and two patient interfaces. The PEP T-piece is supplied with a mouthpiece that can easily be exchanged for a face mask if required. A pressure manometer can be connected to the system using the PEP manometer connector (1956020) designed to be fitted between the T-piece and PEP resistor.



Code	Description	Colour	Box Qty.
The PEP resistors are packaged in bags of 10 units			
8020060	PEP Resistor 6.0mm, purple, single patient use		100
8020050	PEP Resistor 5.0mm, brown, single patient use		100
8020045	PEP Resistor 4.5mm, red, single patient use		100
8020040	PEP Resistor 4.0mm, orange, single patient use		100
8020035	PEP Resistor 3.5mm, green, single patient use		100
8020030	PEP Resistor 3.0mm, blue, single patient use		100
8020025	PEP Resistor 2.5mm, yellow, single patient use		100
8020020	PEP Resistor 2.0mm, white, single patient use		100
8020015	PEP Resistor 1.5mm, black, single patient use		100
PEP T-piece, patient interfaces and connectors			
1956014	PEP T-piece, Inhalation port -22M/15F-22M with one-way valve and mouthpiece, 22F, single patient use		80
1956020	PEP manometer connector, 15M, 22M/15F, 6mm elbow with tube, 1.3m		40
7296001	ClearLite, anaesthetic face mask, size 6, extra large adult, red seal, no hook ring, 22F		20
7295001	ClearLite, anaesthetic face mask, size 5, large adult, orange seal, no hook ring, 22F		35
7294001	ClearLite, anaesthetic face mask, size 4, adult, green seal, no hook ring, 22F		35
7293001	ClearLite, anaesthetic face mask, size 3, small adult, yellow seal, no hook ring, 22F		35
7292001	ClearLite, anaesthetic face mask, size 2, paediatric, white seal, no hook ring, 22F		25
1691040	Nose clip		35

IS10.37 • Issue 6 09.25



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

Please think before you print
Save energy and paper.
If you must print this information sheet please print it double sided.