




# BiTrac™ MaxShield Select™

Total face mask for Non-invasive Ventilation



Critical Care ▪ Pulmonary®



Quality, innovation and choice

Interact with us



[www.intersurgical.com](http://www.intersurgical.com)



## BiTrac™ MaxShield Select™

Our total face mask, the BiTrac MaxShield Select, avoids the nose and disperses pressure over a larger surface area while reducing claustrophobia for the patient. Our interchangeable elbows greatly reduce waste since you can easily remove and replace the mask's elbow when transitioning patients between standard ventilators and bilevel machines.



**3 sizes available**  
from extra large adult  
to small adult

**Breathable headgear**

**Bayonet style lock**  
for secure elbow connection

**Secure clips**

**Interchangeable elbows available**  
to fit your patient's needs

**Soft silicone seal**



### Interchangeable elbows

Code	Description	Box Qty.
313-9599	BiTrac Select interchangeable standard non-vented elbow	10
313-9593	BiTrac Select interchangeable vented elbow with anti-asphyxiation valve	10
313-9598	BiTrac Select interchangeable non-vented elbow with anti-asphyxiation valve	10
313-9596	BiTrac Select interchangeable non-vented elbow with anti-asphyxiation valve and bronchoscopy port	10
313-9597	BiTrac Select interchangeable non-vented elbow with bronchoscopy port	10



### BiTrac™ MaxShield Select™ total face mask with interchangeable vented elbow with anti-asphyxiation valve

Code	Description	Box Qty.
313-9555	BiTrac MaxShield Select total face mask with interchangeable vented elbow with anti-asphyxiation valve, extra large adult	10
313-9552	BiTrac MaxShield Select total face mask with interchangeable vented elbow with anti-asphyxiation valve, large adult	10
313-9554	BiTrac MaxShield Select total face mask with interchangeable vented elbow with anti-asphyxiation valve, small adult	10



### BiTrac MaxShield Select total face mask with interchangeable standard non-vented elbow

Code	Description	Box Qty.
313-9556	BiTrac MaxShield Select total face mask with interchangeable standard non-vented elbow, extra large adult	10
313-9553	BiTrac MaxShield Select total face mask with interchangeable standard non-vented elbow, large adult	10
313-9551	BiTrac MaxShield Select total face mask with interchangeable standard non-vented elbow, small adult	10




### BiTrac MaxShield Select total face mask with interchangeable non-vented elbow with anti-asphyxiation valve

Code	Description	Box Qty.
313-9559	BiTrac MaxShield Select total face mask with interchangeable non-vented elbow with anti-asphyxiation valve, extra large adult	10
313-9558	BiTrac MaxShield Select total face mask with interchangeable non-vented elbow with anti-asphyxiation valve, large adult	10
313-9557	BiTrac MaxShield Select total face mask with interchangeable non-vented elbow with anti-asphyxiation valve, small adult	10




**BiTrac™ MaxShield Select™ paediatric total face mask with interchangeable vented elbow with anti-asphyxiation valve**

Code	Description	Box Qty.
313-9527	BiTrac MaxShield Select paediatric total face mask with interchangeable vented elbow with anti-asphyxiation valve, 15 cm extender tube with 22M swivel connection, extra small	10
 313-9526	BiTrac MaxShield Select paediatric total face mask with interchangeable vented elbow with anti-asphyxiation valve, 15 cm extender tube with 22M swivel connection, extra extra small	10



**BiTrac MaxShield Select paediatric total face mask with interchangeable standard non-vented elbow**

Code	Description	Box Qty.
 313-9525	BiTrac MaxShield Select paediatric total face mask with interchangeable standard non-vented elbow, extra small	10
313-9524	BiTrac MaxShield Select paediatric total face mask with interchangeable standard non-vented elbow, extra extra small	10

**Replacement headgears**

Code	Description	Box Qty.
313-8344	Replacement headgear for BiTrac MaxShield, size extra large and large	10
313-8347	Replacement headgear for BiTrac MaxShield, size small	10
313-9529	Replacement headgear for paediatric BiTrac MaxShield, size extra small	10
313-9528	Replacement headgear for paediatric BiTrac MaxShield, size extra extra small	10

 illustrated

BiTrac MaxShield Select is manufactured by Pulmodyne Inc, USA and distributed by Intersurgical Ltd

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

IS1.60 ▪ Issue 1 05.22



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com